

TRIB2 Antibody, Biotin conjugated



PACO53505

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Interacts with MAPK kinases and regulates activation of MAP kinases. Does not display kinase activity.

Gene ID:

TRIB2

Uniprot

Q92519

Synonyms:

Tribbles homolog 2 (TRB-2), TRIB2, TRB2

Immunogen:

Recombinant Human Tribbles homolog 2 protein (204-337AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A