

# ACY3 Antibody, HRP conjugated



PACO53515

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Plays an important role in deacetylating mercapturic acid, in kidney proximal tubules. Also acts on N-acetyl-aromatic amino acid, .

**Gene ID:**

ACY3

**Uniprot**

Q96HD9

**Synonyms:**

N-acyl-aromatic-L-amino acid, amidohydrolase (carboxylate-forming) (EC 3.5.1.114) (Acylase III) (Aminoacylase-3) (ACY-3) (Aspartoacylase-2) (Hepatitis C virus core-binding protein 1) (HCBP1) (HCV core-binding protein 1), ACY3, ASPA2

**Immunogen:**

Recombinant Human N-acyl-aromatic-L-amino acid, amidohydrolase (carboxylate-forming) protein (5-319AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A