

ACY3 Antibody, Biotin conjugated



PACO53517

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays an important role in deacetylating mercapturic acid, in kidney proximal tubules. Also acts on N-acetyl-aromatic amino acid, .

Gene ID:

ACY3

Uniprot

Q96HD9

Synonyms:

N-acyl-aromatic-L-amino acid, amidohydrolase (carboxylate-forming) (EC 3.5.1.114) (Acylase III) (Aminoacylase-3) (ACY-3) (Aspartoacylase-2) (Hepatitis C virus core-binding protein 1) (HCBP1) (HCV core-binding protein 1), ACY3, ASPA2

Immunogen:

Recombinant Human N-acyl-aromatic-L-amino acid, amidohydrolase (carboxylate-forming) protein (5-319AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A