

SCG3 Antibody



PACO53518

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

cytoplasmic vesicle, extracellular region, secretory granule lumen, poly(A) RNA binding, platelet degranulation.

Gene ID:

SCG3

Uniprot

Q8WXD2

Synonyms:

Secretogranin-3 (Secretogranin III) (SgIII), SCG3


Immunogen:

Recombinant Human Secretogranin-3 protein (205-432AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunofluorescent analysis of HeLa cells using PACO53518 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).