

PREB Antibody, HRP conjugated



PACO53551

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Guanine nucleotide exchange factor that specifically activates the small GTPase SAR1B. Mediates the recruitment of SAR1B and other COPII coat components to endoplasmic reticulum membranes and is therefore required for the formation of COPII transport vesicles from the ER.

Gene ID:

PREB

Uniprot

Q9HCU5

Synonyms:

Prolactin regulatory element-binding protein (Mammalian guanine nucleotide exchange factor mSec12), PREB, SEC12

Immunogen:

Recombinant Human Prolactin regulatory element-binding protein (98-283AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A