

NDUFA4L2 Antibody, HRP conjugated



PACO53559

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

mitochondrial respiratory chain complex IV, NADH dehydrogenase (ubiquinone) activity.

Gene ID:

NDUFA4L2

Uniprot

Q9NRX3

Synonyms:

NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 4-like 2 (NADH-ubiquinone oxidoreductase MLRQ subunit homolog) (NUOMS), NDUFA4L2

Immunogen:

Recombinant Human NADH dehydrogenase [ubiquinone] 1 α subcomplex subunit 4-like 2 protein (48-87AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A