

PRDM13 Antibody



PACO53566

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

PRDM13

Uniprot

Q9H4Q3

Synonyms:

PR domain zinc finger protein 13 (EC 2.1.1) (PR domain-containing protein 13), PRDM13, PFM10


Immunogen:

Recombinant Human PR domain zinc finger protein 13 protein (2-255AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: Mouse heart tissue. All lanes: PRDM13 antibody at 2.7 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 74 kDa. Observed band size: 74 kDa.