

# ACP6 Antibody, FITC conjugated



PACO53572

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Hydrolyzes lysophosphatidic acid, (LPA) containing a medium length fatty acid, chain to the corresponding monoacylglycerol. Has highest activity with lysophosphatidic acid, containing myristate (C14:0), monounsaturated oleate (C18:1) or palmitate (C16:0), and lower activity with C18:0 and C6:0 lysophosphatidic acid, .

**Gene ID:**

ACP6

**Uniprot**

Q9NPH0

**Synonyms:**

Lysophosphatidic acid, phosphatase type 6 (EC 3.1.3.2) (acid, phosphatase 6, lysophosphatidic) (acid, phosphatase-like protein 1) (PACPL1), ACP6, ACPL1 LPAP

**Immunogen:**

Recombinant Human Lysophosphatidic acid, phosphatase type 6 protein (35-324AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A