

BTAf1 Antibody, HRP conjugated



PACO53611

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Regulates transcription in association with TATA binding protein (TBP). Removes TBP from the TATA box in an ATP-dependent manner.

Gene ID:

BTAf1

Uniprot

O14981

Synonyms:

TATA-binding protein-associated factor 172 (EC 3.6.4) (ATP-dependent helicase BTAf1) (B-TFIID transcription factor-associated 170 kDa subunit) (TAF(II)170) (TBP-associated factor 172) (TAF-172), BTAf1, TAF172

Immunogen:

Recombinant Human TATA-binding protein-associated factor 172 protein (1750-1849AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A