

# GPS2 Antibody, FITC conjugated



PACO53668

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Suppresses G-protein- and mitogen-activated protein kinase-mediated signal transduction.

**Gene ID:**

GPS2

**Uniprot**

Q13227

**Synonyms:**

G protein pathway suppressor 2 (GPS-2), GPS2

**Immunogen:**

Recombinant Human G protein pathway suppressor 2 protein (1-327AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A