

# tal1 Antibody, FITC conjugated



PACO53780

---

## Product Information

**Size:**

50ug

**Reactivity:**

Zebrafish

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

tal1

**Uniprot**

O93507

**Synonyms:**

T-cell acute lymphocytic leukemia protein 1 homolog (TAL-1) (Stem cell protein) (zSCL), tal1, scl tal-1

**Immunogen:**

Recombinant Zebrafish T-cell acute lymphocytic leukemia protein 1 homolog protein (1-324AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A