

PACO53802

---

## Product Information

**Size:**

50ug

**Reactivity:**

Bacteroides fragilis

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:5000

**Protein Background:**

Diarrheal toxin that hydrolyzes gelatin, azocoll, actin, tropomyosin, and fibrinogen.

**Gene ID:**

btfP

**Uniprot**

P54355

**Synonyms:**

Fragilysin (EC 3.4.24.74) (Enterotoxin), btFP

**Immunogen:**


Recombinant Bacteroides fragilis Fragilysin protein (26-405AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Western Blot. Positive WB detected in Recombinant protein. All lanes: btfP antibody at 3 $\mu$ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 58 kDa. Observed band size: 58 kDa.