

THI2.1 Antibody



PACO53810

Product Information

Size:

50ug

Reactivity:

Viscum album

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Thionins are small plant proteins which are toxic to animal cells. They seem to exert their toxic effect at the level of the cell membrane. Their precise function is not known.

Gene ID:

THI2.1

Uniprot

P01538

Synonyms:

Viscotoxin-A3 [Cleaved into: Viscotoxin-A3; acid, c protein], THI2.1


Immunogen:

Recombinant Viscum album Viscotoxin-A3 protein (27-72AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: TH12.1 antibody at 2.7 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 20 kDa. Observed band size: 20 kDa.