

# TSR1 Antibody, FITC conjugated



PACO53860

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Required during maturation of the 40S ribosomal subunit in the nucleolus.

**Gene ID:**

TSR1

**Uniprot**

Q2NL82

**Synonyms:**

Pre-rRNA-processing protein TSR1 homolog, TSR1, KIAA1401

**Immunogen:**

Recombinant Human Pre-rRNA-processing protein TSR1 homolog protein (455-804AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A