

VPS16 Antibody, FITC conjugated



PACO53900

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	VPS16
Source:	Uniprot
Rabbit	Q6GMQ7
Isotype:	Synonyms:
IgG	Vacuolar protein sorting-associated protein 16 homolog (hVPS16), VPS16
Applications:	Immunogen:
ELISA	Recombinant Human VPS16 protein (7-89AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A