

# MKRN2 Antibody, Biotin conjugated



PACO53921

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

E3 ubiquitin ligase catalyzing the covalent attachment of ubiquitin moieties onto substrate proteins.

**Gene ID:**

MKRN2

**Uniprot**

Q9H000

**Synonyms:**

Probable E3 ubiquitin-protein ligase makorin-2 (EC 2.3.2.27) (RING finger protein 62) (RING-type E3 ubiquitin transferase makorin-2), MKRN2, RNF62

**Immunogen:**

Recombinant Human Probable E3 ubiquitin-protein ligase makorin-2 protein (29-170AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A