

YBEY Antibody



PACO53970

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

nucleoplasm, nucleus.

Gene ID:

YBEY

Uniprot

P58557

Synonyms:

Endoribonuclease YbeY (EC 3.1), YBEY, C21orf57


Immunogen:

Recombinant Human Endoribonuclease YbeY protein (54-167AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry analysis of human colon cancer using PACO53970 at dilution of 1:100.



Immunocytochemistry analysis of human tonsil tissue using PACO53970 at dilution of 1:100.