

# ABR Antibody, FITC conjugated



PACO54120

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

GTPase-activating protein for RAC and CDC42. Promotes the exchange of RAC or CDC42-bound GDP by GTP, thereby activating them.

**Gene ID:**

ABR

**Uniprot**

Q12979

**Synonyms:**

Active breakpoint cluster region-related protein, ABR

**Immunogen:**

Recombinant Human Active breakpoint cluster region-related protein (797-859AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A