

# HLA-B Antibody, HRP conjugated



PACO54123

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Involved in the presentation of foreign antigens to the immune system.

**Gene ID:**

HLA-B

**Uniprot**

Q29718

**Synonyms:**

HLA class I histocompatibility antigen, B-82 alpha chain (MHC class I antigen B\*82), HLA-B, HLAB

**Immunogen:**

Recombinant Human HLA class I histocompatibility antigen, B-82 &alpha; chain protein (197-306AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A