

# WWC2 Antibody, HRP conjugated



PACO54131

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

WWC2

**Uniprot**

Q6AWC2

**Synonyms:**

Protein WWC2 (BH-3-only member B) (WW domain-containing protein 2), WWC2, BOMB

**Immunogen:**

Recombinant Human Protein WWC2 protein (828-1127AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A