

RAD9A Antibody



PACO54178

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

RAD9A

Uniprot

Q99638

Synonyms:

Cell cycle checkpoint control protein RAD9A (hRAD9) (EC 3.1.11.2) (DNA repair exonuclease rad9 homolog A), RAD9A


Immunogen:

Recombinant Human Cell cycle checkpoint control protein RAD9A protein (267-341AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry analysis of human kidney tissue using PACO54178 at dilution of 1:100.



Immunocytochemistry analysis of human breast cancer using PACO54178 at dilution of 1:100.