

RAD9A Antibody, FITC conjugated



PACO54180

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

RAD9A

Uniprot

Q99638

Synonyms:

Cell cycle checkpoint control protein RAD9A (hRAD9) (EC 3.1.11.2) (DNA repair exonuclease rad9 homolog A), RAD9A

Immunogen:

Recombinant Human Cell cycle checkpoint control protein RAD9A protein (267-341AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A