

# PSORS1C1 Antibody, HRP conjugated



PACO54215

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:****Gene ID:**

PSORS1C1

**Uniprot**

Q9UIG5

**Synonyms:**

Psoriasis susceptibility 1 candidate gene 1 protein (Protein SEEK1), PSORS1C1, C6orf16  
SEEK1

**Immunogen:**

Recombinant Human Psoriasis susceptibility 1 candidate gene 1 protein (1-152AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A