

PACO54250

Product Information

Size:

50ug

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Required for pectoral fin formation. Together with tbx5b, involved in eye and heart development. Required for the looping stage of heart development.

Gene ID:

tbx5a

Uniprot

Q9IAK8

Synonyms:

T-box transcription factor TBX5-A (zftbx5a) (T-box protein 5) (zTbx5), tbx5a, heartstrings
hst tbx5 tbx5.1

Immunogen:

Recombinant Zebrafish T-box transcription factor TBX5-A protein (1-492AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A