

FOSB Antibody, Biotin conjugated



PACO54297

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

FosB interacts with Jun proteins enhancing their DNA binding activity.

Gene ID:

FOSB

Uniprot

P53539

Synonyms:

Protein fosB (G0/G1 switch regulatory protein 3), FOSB, G0S3

Immunogen:

Recombinant Human Protein fosB protein (1-41AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A