

GRB10 Antibody, FITC conjugated



PACO54304

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

GRB10

Uniprot

Q13322

Synonyms:

Growth factor receptor-bound protein 10 (GRB10 adapter protein) (Insulin receptor-binding protein Grb-IR), GRB10, GRBIR KIAA0207

Immunogen:

Recombinant Human Growth factor receptor-bound protein 10 protein (15-180AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A