

IPP Antibody, FITC conjugated



PACO54312

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play a role in organizing the actin cytoskeleton.

Gene ID:

IPP

Uniprot

Q9Y573

Synonyms:

Actin-binding protein IPP (Intracisternal A particle-promoted polypeptide) (IPP) (Kelch-like protein 27), IPP, KLHL27

Immunogen:

Recombinant Human Actin-binding protein IPP protein (444-515AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A