

LCT Antibody, Biotin conjugated



PACO54325

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

LPH splits lactose in the small intestine.

Gene ID:

LCT

Uniprot

P09848

Synonyms:

Lactase-phlorizin hydrolase (Lactase-glycosylceramidase) [Includes: Lactase (EC 3.2.1.108); Phlorizin hydrolase (EC 3.2.1.62)], LCT, LPH

Immunogen:

Recombinant Human Lactase-phlorizin hydrolase protein (983-1252AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A