

NOTCH1 Antibody, FITC conjugated



PACO54348

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

NOTCH1

Uniprot

P46531

Synonyms:

Neurogenic locus notch homolog protein 1 (Notch 1) (hN1) (Translocation-associated notch protein TAN-1) [Cleaved into: Notch 1 extracellular truncation (NEXT); Notch 1 intracellular domain (NICD)], NOTCH1, TAN1

Immunogen:

Recombinant Human Neurogenic locus notch homolog protein 1 protein (2428-2555AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A