

# OLIG1 Antibody, Biotin conjugated



PACO54353

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Promotes formation and maturation of oligodendrocytes, especially within the brain. Cooperates with OLIG2 to establish the pMN domain of the embryonic neural tube.

**Gene ID:**

OLIG1

**Uniprot**

Q8TAK6

**Synonyms:**

Oligodendrocyte transcription factor 1 (Oligo1) (Class B basic helix-loop-helix protein 6) (bHLHb6) (Class E basic helix-loop-helix protein 21) (bHLHe21), OLIG1, BHLHB6 BHLHE21

**Immunogen:**

Recombinant Human Oligodendrocyte transcription factor 1 protein (17-105AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A