

OLIG1 Antibody, HRP conjugated



PACO54355

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Promotes formation and maturation of oligodendrocytes, especially within the brain. Cooperates with OLIG2 to establish the pMN domain of the embryonic neural tube.

Gene ID:

OLIG1

Uniprot

Q8TAK6

Synonyms:

Oligodendrocyte transcription factor 1 (Oligo1) (Class B basic helix-loop-helix protein 6) (bHLHb6) (Class E basic helix-loop-helix protein 21) (bHLHe21), OLIG1, BHLHB6 BHLHE21

Immunogen:

Recombinant Human Oligodendrocyte transcription factor 1 protein (17-105AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A