

RBP1 Antibody



PACO54382

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Intracellular transport of retinol.

Gene ID:

RBP1

Uniprot

P09455

Synonyms:

Retinol-binding protein 1 (Cellular retinol-binding protein) (CRBP) (Cellular retinol-binding protein I) (CRBP-I), RBP1, CRBP1


Immunogen:

Recombinant Human Retinol-binding protein 1 protein (70-112AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO54382 at dilution of 1:100.