

# Invasin Antibody



PACO54458

---

## Product Information

**Size:**

50ug

**Reactivity:**

Yersinia enterocolitica

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Invasin is a protein that allows enteric bacteria to penetrate cultured mammalian cells. The entry of invasin in the cell is mediated by binding several beta-1 chain integrins.

**Gene ID:****Uniprot**

P19196

**Synonyms:**

Invasin

**Immunogen:**

Recombinant Yersinia enterocolitica Invasin protein (651-835AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A