

Invasin Antibody, Biotin conjugated



PACO54461

Product Information

Size:

50ug

Reactivity:

Yersinia enterocolitica

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Invasin is a protein that allows enteric bacteria to penetrate cultured mammalian cells. The entry of invasin in the cell is mediated by binding several beta-1 chain integrins.

Gene ID:**Uniprot**

P19196

Synonyms:

Invasin

Immunogen:

Recombinant Yersinia enterocolitica Invasin protein (651-835AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A