

TBX2 Antibody, Biotin conjugated



PACO54473

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in the transcriptional regulation of genes required for mesoderm differentiation. Probably plays a role in limb pattern formation. Acts as a negative regulator of PML function in cellular senescence. May be required for cardiac atrioventricular canal formation.

Gene ID:

TBX2

Uniprot

Q13207

Synonyms:

T-box transcription factor TBX2 (T-box protein 2), TBX2

Immunogen:

Recombinant Human T-box transcription factor TBX2 protein (299-573AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A