

PKD2 Antibody, Biotin conjugated



PACO54485

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PKD2

Uniprot

Q13563

Synonyms:

Polycystin-2 (PC2) (Autosomal dominant polycystic kidney disease type II protein) (Polycystic kidney disease 2 protein) (Polycystin) (R48321) (Transient receptor potential cation channel subfamily P member 2), PKD2, TRPP2

Immunogen:

Recombinant Human Polycystin-2 protein (680-968AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A