

PAG 1 Antibody, FITC conjugated



PACO54500

Product Information

Size:

50ug

Reactivity:

Bovine

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Appears to be proteolytically inactive.

Gene ID:**Uniprot**

Q29432

Synonyms:

Pregnancy-associated glycoprotein 1 (PAG 1) (EC 3.4.23) (Pregnancy-specific protein B) (PSP-B)

Immunogen:

Recombinant Bovine Pregnancy-associated glycoprotein 1 protein (54-380AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A