

# TOP2B Antibody, Biotin conjugated



PACO54517

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Control of topological states of DNA by transient breakage and subsequent rejoining of DNA strands. Topoisomerase II makes double-strand breaks.

**Gene ID:**

TOP2B

**Uniprot**

Q02880

**Synonyms:**

DNA topoisomerase 2-beta (EC 5.99.1.3) (DNA topoisomerase II, beta isozyme), TOP2B

**Immunogen:**

Recombinant Human DNA topoisomerase 2-&beta; protein (1400-1626AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A