

# ADGRA2 Antibody, FITC conjugated



PACO54556

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Endothelial receptor which functions as a WNT7-specific coactivator of canonical Wnt signaling. Required for normal endothelial cell sprouting and migration in the forebrain and neural tube. Has a major role in establishing the blood-brain barrier. Binds to the glycosaminoglycans heparin, heparin sulfate, chondroitin sulfate and dermatan sulfate.

**Gene ID:**

ADGRA2

**Uniprot**

Q96PE1

**Synonyms:**

Adhesion G protein-coupled receptor A2 (G-protein coupled receptor 124) (Tumor endothelial marker 5), ADGRA2, GPR124 KIAA1531 TEM5

**Immunogen:**

Recombinant Human Adhesion G protein-coupled receptor A2 protein (1069-1338AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A