

RBPJL Antibody, HRP conjugated



PACO54575

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Putative transcription factor, which cooperates with EBNA2 to activate transcription.

Gene ID:

RBPJL

Uniprot

Q9UBG7

Synonyms:

Recombining binding protein suppressor of hairless-like protein (Transcription factor RBP-L), RBPJL, RBPL, RBPSUHL

Immunogen:

Recombinant Human Recombining binding protein suppressor of hairless-like protein (1-151AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A