

ABCB4 Antibody



PACO54602

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:200-1:500

Protein Background:

Gene ID:

ABCB4

Uniprot

P21439

Synonyms:

Phosphatidylcholine translocator ABCB4 (ATP-binding cassette sub-family B member 4) (Multidrug resistance protein 3) (EC 3.6.3.44) (P-glycoprotein 3), ABCB4, MDR3 PGY3


Immunogen:

Recombinant Human Phosphatidylcholine translocator ABCB4 protein (530-693AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human melanoma using PACO54602 at dilution of 1:100.



Immunofluorescence staining of A549 cells with PACO54602 at 1:200, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).