

ABCB7 Antibody, HRP conjugated



PACO54607

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Could be involved in the transport of heme from the mitochondria to the cytosol. Plays a central role in the maturation of cytosolic iron-sulfur (Fe/S) cluster-containing proteins.

Gene ID:

ABCB7

Uniprot

O75027

Synonyms:

ATP-binding cassette sub-family B member 7, mitochondrial (ATP-binding cassette transporter 7) (ABC transporter 7 protein), ABCB7, ABC7

Immunogen:

Recombinant Human ATP-binding cassette sub-family B member 7, mitochondrial protein (23-140AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A