

ABCD1 Antibody



PACO54614

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Probable transporter. The nucleotide-binding fold acts as an ATP-binding subunit with ATPase activity.

Gene ID:

ABCD1

Uniprot

P33897

Synonyms:

ATP-binding cassette sub-family D member 1 (Adrenoleukodystrophy protein) (ALDP), ABCD1, ALD


Immunogen:

Recombinant Human ATP-binding cassette sub-family D member 1 protein (632-745AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human lung cancer using PACO54614 at dilution of 1:100.