

# ACADVL Antibody, FITC conjugated



PACO54620

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Active toward esters of long-chain and very long chain fatty acid, such as palmitoyl-CoA, mysritoyl-CoA and stearoyl-CoA. Can accommodate substrate acyl chain lengths as long as 24 carbons, but shows little activity for substrates of less than 12 carbons.

**Gene ID:**

ACADVL

**Uniprot**

P49748

**Synonyms:**

Very long-chain specific acyl-CoA dehydrogenase, mitochondrial (VLCAD) (EC 1.3.8.9), ACADVL, VLCAD

**Immunogen:**

Recombinant Human Very long-chain specific acyl-CoA dehydrogenase, mitochondrial protein (403-497AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A