

ACADVL Antibody, Biotin conjugated



PACO54621

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Active toward esters of long-chain and very long chain fatty acid, such as palmitoyl-CoA, mysritoyl-CoA and stearoyl-CoA. Can accommodate substrate acyl chain lengths as long as 24 carbons, but shows little activity for substrates of less than 12 carbons.

Gene ID:

ACADVL

Uniprot

P49748

Synonyms:

Very long-chain specific acyl-CoA dehydrogenase, mitochondrial (VLCAD) (EC 1.3.8.9), ACADVL, VLCAD

Immunogen:

Recombinant Human Very long-chain specific acyl-CoA dehydrogenase, mitochondrial protein (403-497AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A