

ADCY6 Antibody, Biotin conjugated



PACO54629

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

ADCY6

Uniprot

O43306

Synonyms:

Adenylate cyclase type 6 (EC 4.6.1.1) (ATP pyrophosphate-lyase 6) (Adenylate cyclase type VI) (Adenylyl cyclase 6) (Ca(2+)-inhibitabile adenylyl cyclase), ADCY6, KIAA0422

Immunogen:

Recombinant Human Adenylate cyclase type 6 protein (9-118AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A