

CYP2B6 Antibody, Biotin conjugated



PACO54733

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Cytochromes P450 are a group of heme-thiolate monooxygenases. In liver microsomes, this enzyme is involved in an NADPH-dependent electron transport pathway. It oxidizes a variety of structurally unrelated compounds, including steroids, fatty acid, , and xenobiotics. Acts as a 1,4-cineole 2-exo-monooxygenase.

Gene ID:

CYP2B6

Uniprot

P20813

Synonyms:

Cytochrome P450 2B6 (EC 1.14.13) (1,4-cineole 2-exo-monooxygenase) (CYP11B6) (Cytochrome P450 I1B1), CYP2B6

Immunogen:

Recombinant Human Cytochrome P450 2B6 protein (250-343AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A