

PAC054734

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Catalyzes the conversion of dihydroorotate to orotate with quinone as electron acceptor.

Gene ID:

DHODH

Uniprot

Q02127

Synonyms:

Dihydroorotate dehydrogenase (quinone), mitochondrial (DHOdehase) (EC 1.3.5.2) (Dihydroorotate oxidase), DHODH


Immunogen:

Recombinant Human Dihydroorotate dehydrogenase (quinone), mitochondrial protein (241-395AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO54734 at dilution of 1:100.