

MAST3 Antibody, FITC conjugated



PACO54872

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

intracellular, protein serine/threonine kinase activity, cytoskeleton organization, intracellular signal transduction, peptidyl-serine phosphorylation.

Gene ID:

MAST3

Uniprot

O60307

Synonyms:

Microtubule-associated serine/threonine-protein kinase 3 (EC 2.7.11.1), MAST3, KIAA0561

Immunogen:

Recombinant Human Microtubule-associated serine/threonine-protein kinase 3 protein (1101-1309AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A