

NCAM2 Antibody, Biotin conjugated



PACO54893

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play important roles in selective fasciculation and zone-to-zone projection of the primary olfactory axons.

Gene ID:

NCAM2

Uniprot

O15394

Synonyms:

Neural cell adhesion molecule 2 (N-CAM-2) (NCAM-2), NCAM2, NCAM21

Immunogen:

Recombinant Human Neural cell adhesion molecule 2 protein (209-333AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A